

# Medical Protective Mask (Ear-hanging style)

## 一、Product name, model and specification

Product name: Medical Protective Mask

Model: folding type (Non-sterile)、folding type (sterile);

Specification: L 16cm×10cm、M 16cm×9cm、S 14cm×9cm

Model: arch form (Non-sterile)、arch form (sterile);

Specification: L 14cm×11cm、M 13cm×13cm、S 13cm×12cm

## 二、Name, address, production address, contact information and after-sales service unit of the registered person/manufacturer

- 1、Name: Hunan Zhenheyikang Medical Instrument Co.,LTD
- 2、Address: No.6 Building Jingxiang Energy, No.55 Xiaguang east road, Gaoxin district Xiangtan, Hunan
- 3、Production address: NO.6 Building, NO.8 building, 3<sup>rd</sup> and 4<sup>th</sup> floor Jingxiang Energy, No.55 Xiaguang east road, Gaoxin district Xiangtan, Hunan, NO 19 Chutian Road, Yuetang District, Xiangtan City, Hunan
- 4、Contact number: +86 17711627084
- 5、After-sales service unit: Hunan Zhenheyikang Medical Instrument Co.,LTD

三、Manufacturing License NO: Hunan Medical Products Administration Production Licence No. 20150013

四、Product Registration Certificate No, Product Technical Requirements No: HNFDA Certified 20202140323

## 五、Structure and composition, scope of application

It is composed of mask body, nose clip and mask belt. The inner body and outer layers of the mask are ordinary non-woven cloth, and the middle layer is melt-blown cloth. The nose clips are made of polypropylene material, and the mask bands are made of elastic bands. Non-sterile supply.

Suitable for hygiene care in a general medical setting without body fluids and splash risk.

六、Use method: open the bag can be used

## 七、Contraindications, precautions, warnings and reminders

Product contraindications: None

Damaged packaging prohibited use, expired products prohibited use;

Disposable use, please destroy after use

Avoid high temperature, humidity, direct sunlight

## 八、Interpretation of graphics, symbols, abbreviations and other contents used in medical device labels



一次性使用



参考使用说明书



包装破损切勿使用



未灭菌

## 九、Installation and use instructions or illustrations

It is recommended to use no more than 4 hours.

Instructions for use: 1. Open the package; 2. Hold the end of the nose clip and hold the edge of the ear strap with both hands to open the mask; 3. Hang the ear strap on the ear. Consider the difference in faces, ear hanging with plastic hooks. Adjust and wear, refer to the figure below.



## 十、Product maintenance and maintenance methods, special storage, transportation conditions, methods

Products should be stored in relative humidity  $\leq 80\%$ , non-corrosive gas and well-ventilated indoor, ordinary transport, in the transportation, storage process should be kept clean and pollution-free, strictly prohibited rain and snow immersion.

## 十一、Date of production, shelf life or expiration date

Valid for three years, production batch number, production date, expiry date please refer to the small package or qualification label.

**十二、 Accessories list, including accessories, appendix, lossitage replacement cycle and replacement method description**

This product has n o accessories

**十三、 Date and version of preparation or revision of the specification**

This specification was revised in February 2023, version 2023.

。