

Non-woven and Seam technology

Non-woven at both sides provides cloth like feeling to human skin. Micro-porous film provides breathable but liquid&blood resistance property. Conforms to EN13795-High performance (Critical Areal) and is recommended to use for surgical gowns and drapes.

Features

- *Conforms to EN13795-High performance(Critical Areal)
- *High MVTR(over 4000gsm.24hrs)
- *ETO sterilization available
- *ASTM F1670 passed for blood resistance.



Stitched Seam

Added strength at the seams and more protection against particles



Ultrasonic Seam

Cost effective seams with ideal protection against particles and liquids



Bound Stitched Seam

More strength at the seams and higher protection against particles and liquids



Stitched with Tape Seam

Ultimate protection against particles and liquids.