

Medical Elastic Bandage

一、Product name, model and specification

Product name : Medical Elastic Bandage

Model : Abdominal type、Fetal supervision and fixation type、First aid dressing type、conventional type、Neck type、Carpal joint type、Knee type. Specification : Abdominal type、Neck type、Carpal joint type、Knee type : M、L ; Fetal supervision and fixation type、First aid dressing type、

Conventional type : 1~110cm×100~1000cm ; Neck type、Carpal joint type、Knee type : S 、 M 、 L.

二、Name、address、production address、contact information and after-sales service unit of the registered person/manufacturer

1、Name: Hunan Zhenheyikang Medical Instrument Co.,LTD

2、Address: No.6 Building Jingxiang Energy, No.55 Xianguang east road, Gaoxin district Xiangtan, Hunan

3、Production address:NO.6 Building, NO.8 Building,5th floor Jingxiang Energy, No.55 Xianguang east road,Gaoxin district Xiangtan, Hunan

4、Contact number: +86 17711627084

5、After-sales service unit: Hunan Zhenheyikang Medical Instrument Co.,LTD

三、Production Registration Certificate No: Xiangtan Medical Products Administration Production Registration NO 20140004

四、Product Registration Certificate No、Product Technical Requirements No: XTFDA Certified 20190012

五、Product description, intended use

Medical elastic bandage is made of spandex, polyester, human cotton, pure cotton, mesh cloth, buttons, self-adhesive cloth, water absorption pad, medical skim cotton, medical gauze sheet, plastic film, chemical fiber and other materials, through weaving, elastic, sewing and refined components. The product has elastic or self-adhesive characteristics. The product is banded or tubular, elastic and without direct contact with the wound surface. Non-sterile provided, single use.

It is used to provide binding force to the wound dressing or limb to play the role of dressing and fixing.

六、Use method: Open the bag can be used, directly to treat the wound

七、Contraindications, precautions, warnings and reminders

Product contraindications: None

Damaged packaging prohibited use, expired products prohibited use

Disposable use, please destroy after use

Avoid high temperature, humidity, direct sunlight

八、Interpretation of graphics, symbols, abbreviations and other contents used in medical device labels



一次性使用



参考使用说明书



包装破损切勿使用

九、Installation and use instructions or illustrations

This product does not need to be installed, open the bag can be used.

十、Product maintenance and maintenance methods, special storage, transportation conditions, methods

Products should be stored in a non-corrosive gas and well-ventilated room, transported by ordinary transportation means, keep clean and pollution-free, rain and snow immersion.

十一、Date of production, shelf life or expiration date

Valid for four years, production batch number, production date, expiry date please refer to the small package or qualification label.

十二、Accessories list, including accessories, appendix, lossitage replacement cycle and replacement method description

This product has no accessories

十三、Date and version of preparation or revision of the specification

This specification was revised in March 2023, version 2023.