

Medical Protective Mask (Ear-hanging style)

一、Product name, model and specification

Product name: Medical Protective Mask

Model: folding type (Non-sterile)、folding type (sterile);

Specification: L 16cm×10cm、M 16cm×9cm、S 14cm×9cm

Model: arch form (Non-sterile)、arch form (sterile);

Specification: L 14cm×11cm、M 13cm×13cm、S 13cm×12cm

二、Name, address, production address, contact information and after-sales service unit of the registered person/manufacturer

- 1、Name: Hunan Zhenheyikang Medical Instrument Co.,LTD
- 2、Address: No.6 Building Jingxiang Energy, No.55 Xiaguang east road, Gaoxin district Xiangtan, Hunan
- 3、Production address: NO.6 Building, NO.8 building, 3rd and 4th floor Jingxiang Energy, No.55 Xiaguang east road, Gaoxin district Xiangtan, Hunan, NO 19 Chutian Road, Yuetang District, Xiangtan City, Hunan
- 4、Contact number: +86 17711627084
- 5、After-sales service unit: Hunan Zhenheyikang Medical Instrument Co.,LTD

三、Manufacturing License NO: Hunan Medical Products Administration Production Licence No. 20150013

四、Product Registration Certificate No, Product Technical Requirements No: HNFDA Certified 20202140323

五、Structure and composition, scope of application

It is composed of mask body, nose clip and mask belt. The inner body and outer layers of the mask are ordinary non-woven cloth, and the middle layer is melt-blown cloth. The nose clips are made of polypropylene material, and the mask bands are made of elastic bands. Non-sterile supply.

Suitable for hygiene care in a general medical setting without body fluids and splash risk.

六、Use method: open the bag can be used

七、Contraindications, precautions, warnings and reminders

Product contraindications: None

Damaged packaging prohibited use, expired products prohibited use;

Disposable use, please destroy after use

Avoid high temperature, humidity, direct sunlight

八、Interpretation of graphics, symbols, abbreviations and other contents used in medical device labels



一次性使用



参考使用说明书



包装破损切勿使用



未灭菌

九、Installation and use instructions or illustrations

It is recommended to use no more than 4 hours.

Instructions for use: 1. Open the package; 2. Hold the end of the nose clip and hold the edge of the ear strap with both hands to open the mask; 3. Hang the ear strap on the ear. Consider the difference in faces, ear hanging with plastic hooks. Adjust and wear, refer to the figure below.



十、Product maintenance and maintenance methods, special storage, transportation conditions, methods

Products should be stored in relative humidity $\leq 80\%$, non-corrosive gas and well-ventilated indoor, ordinary transport, in the transportation, storage process should be kept clean and pollution-free, strictly prohibited rain and snow immersion.

十一、Date of production, shelf life or expiration date

Valid for three years, production batch number, production date, expiry date please refer to the small package or qualification label.

十二、 Accessories list, including accessories, appendix, lossitage replacement cycle and replacement method description

This product has n o accessories

十三、 Date and version of preparation or revision of the specification

This specification was revised in February 2023, version 2023.

。